#### 1.1.1 Room

Temperature Control	Minimum 5°C (41°F), Maximum 40°C (104°F)
Humidity	20% to 70%
Condition	Good Ventilation
Altitude	Up to 2000 meters
Additional	Access to exhaust system for glovebox vent tubes

#### 1.1.2 Floor

Condition	Flat and level, Able to support equipment weight at least 700 Kg (1600
	lbs.)

# 1.2 Gas Supply

# 1.2.1 Working Gas

Use	Required for operating certain system valves and to create an inert atmosphere inside the glovebox
Gas Types	Helium (He), Nitrogen (N <sub>2</sub> ), or Argon (Ar), gas cylinders (or boil-off from cryogenic liquid)
Quality	99.995% to ultra-high purity, 99.999%
Quantity	A type K or T* cylinder at approximate 2,400 ~ 3,000 psi (165 ~ 205 bar). K and T are supplier-dependent designations, but both are 250 cubic feet of gas.
Regulators	CGA350, Two-stage, capable of a delivery pressure of at least 80psi (~5 bar)
Connection	8 mm compression fitting (Provided by Vigor) to 1/4 "or 3/8 "NPT

#### 1.2.2 Regeneration Gas

Use	Regenerates the purifier by reducing copper oxide back to copper metal.
Gas Types	Inert gas to match your working gas plus 5% Hydrogen, for example; 5% Hydrogen (H <sub>2</sub> ) in Helium (He); 5% Hydrogen (H <sub>2</sub> ) in Nitrogen (N <sub>2</sub> ) 5% Hydrogen (H <sub>2</sub> ) in Argon (Ar)
Qualities	99.995% to ultra-high purity, 99.999%
Quantity	Requires approximately ½ to 1 full type K or T cylinder at 2,400 ~ 3,000 psi (165 ~ 205 bar) for each regeneration cycle
Regulators	CGA580 or CGA350, Depending on the gas supplier, two-stage, capable of a delivery pressure of at least 40 psi (~3 bar)
Connection	8 mm compression fitting (Provided by Vigor) to 1/4 "or 3/8 "NPT

# 1.3 Power Supply (Standard USA Configuration, 220VAC units are also available)

Main Power	Dedicated 110V – 120VAC, 60Hz, 15A Circuit, Grounded (USA) Connection to the supply is via a standard grounded AC power plug.
Heated Antechamber (option)	~120VAC, 60Hz Units with flat metal heaters are 300W/plate. Units have two (FLC 5A) or three plates (FLC 7.5A). Units with IR heating tubes require a maximum FLC of 8.4A. The heater unit connects to the supply via a standard grounded AC power plug.
Freezer (option)	~120VAC, 60Hz, 300W (Maximum FLC 2.5 A) The freezer unit connects to the supply via a standard grounded AC power plug.